

# ECG & Temperature Sensor with BLE Lab

BME554L - Fall 2025 - Palmeri

Dr. Mark Palmeri, M.D., Ph.D.

2025-06-24

## Table of contents

Git Version Control . . . . .	1
Best Coding Practices . . . . .	2
Firmware Functional Specifications . . . . .	2
BLE Server (Mobile App) . . . . .	3
State Diagram . . . . .	4
Testing & Verification . . . . .	4
How to Generate an ECG Signal . . . . .	4
Grading . . . . .	4
How to Ask for Help . . . . .	5
What to Submit . . . . .	5
Resources . . . . .	5
Heart Rate Service (GATT) . . . . .	5

- Fork this repository to your userspace.
- Add Dr. Palmeri as a **Maintainer**.
- Questions should be asked exclusively through GitLab Issues.

## Git Version Control

- Use best practices for version control (branching, commit messages, etc.).
- Do all development on a dedicated branch that is merged into **main** once it is functional.
- Commits should be very specific to the changes/additions you are making to your code. This will help you and others understand what you did and why you did it.
- On a given development branch, try to implement one small piece of functionality at a time, commit it, and then move on to the next piece of functionality.

### ! Important

You do not want one, monolithic git commit right before you submit your project.

## Best Coding Practices

- Use best coding practices throughout the development of your firmware.
- Functions should be short and do one thing. They should return an exit code that is checked in the calling function, indicating success or failure.
- **MACROS!** Avoid hard-coded values in your code.
- Use structs to organize related data.
- Use libraries for code that is self-contained.
- Use the **LOGGING** module to log errors, warnings, information and debug messages.
- You should not have any compiler/build warnings. The CI script will build against v2.9.0 of the Zephyr SDK.

## Firmware Functional Specifications

- Write all firmware using the state machine framework.
  - Do all device initialization in an **INIT** state.
  - Have an **IDLE** state when the device isn't making any measurements.
  - Have an **ERROR** state if any error exit codes are returned from any functions.
    - \* All 4 LEDs should blink at a 50% duty cycle (**ON:OFF** time), in-phase with each other, in the **ERROR** state.
    - \* An error condition should post an error-related event that causes the device to enter the **ERROR** state.
    - \* The error code should specify the error condition that caused the device to enter the **ERROR** state. For example, you may choose to have a bit array that can capture multiple error conditions.
    - \* A BLE notification should be sent with the error code (see BLE custom service/characteristic below).
  - Implement states of your choosing for the following measurements, calculations and BLE communications.
- Have a heartbeat **LED0** that blinks every 1 second with a 50% duty cycle (**ON:OFF** time) in all states.
- Implement functionality to measure a battery voltage (0-3.0 V) using **AIN0**:
  1. When the device first powers on, and then
  2. Every 1 minute thereafter, but only when in the **IDLE** state.

3. You won't actually be connecting a battery to your device; you can use a power support or another voltage source to input a voltage to `AIN0` to simulate a battery level.
- Have the brightness of `LED1` linearly modulated by the percentage of the battery level.
  - Implement functionality to make two measurements after pressing `BUTTON1`:
    1. Read temperature with your [MCP9808](#) sensor (in degrees Celsius).
    2. Calculate the average heart rate (40-200 BPM) using 25-30 seconds of an ECG signal (ranging from -500 - 500 mV, note this is bipolar) from the function generator (see video on how to setup the function generator to output an ECG signal).
  - Pressing `BUTTON1` during the measurements should post an error and go to the `ERROR` state.
  - Blink `LED2` with a 25% duty cycle (`ON:OFF` time) at the average heart rate after the measurements are complete.
  - Have Bluetooth notifications after the measurements are complete and data have been processed, using the BLE services and characteristics described below.
    - Configure the **DIS (Device Information Service)** to report the device model as your Team Name (come up with something fun).
    - Set the **BAS (Battery Service)** to report the battery level of your device. (This isn't actually a battery level, but we're using the `AIN0` measurement as a surrogate for a battery level.)
    - Set the **Heart Rate Service** to report the average heart rate. (See Resources section below.)
    - Setup a custom service with the following custom characteristics:
      - \* **Temperature** for the I2C temperature sensor data in degrees Celcius.
      - \* **Error Code** for the error code that caused the device to enter the `ERROR` state.
  - `BUTTON2` should clear (turn off) a blinking `LED2`, and if `LED2` is not blinking because a measurement hasn't been taken, then it should log a warning (`LOG_WRN()`) as to why it appears nothing happened.
  - `BUTTON3` should be used to reset the device from the `ERROR` state and return to the `IDLE` state.
  - Use timers, kernel events, work queues, threads and any other Zephyr RTOS features as needed to implement the above functionality.

## BLE Server (Mobile App)

- Your device can connect via BLE to a mobile app called [nRF Connect](#).
- This app can be used to read the services and characteristics that your device is advertising.

## State Diagram

- Generate a detailed state diagram that all states, events and actions for your firmware.
- A starter diagram is provided in the `state_diagram.puml` file, along with its rendering below.
- You should add states, events and actions as needed to fully describe the functionality of your firmware.




Figure 1: State Diagram

## Testing & Verification

Complete the testing analysis described in [testing/final\\_project.ipynb](#) to verify the accuracy of your firmware.

## How to Generate an ECG Signal

- [WaveStation 2012 AWG](#)
- [Digilent Waveforms Script](#)

## Grading

- This final project is worth 75% of your grade. *Absolutely no late submissions will be accepted.*
- Git version control will be graded based on best practices.
- Firmware will be graded based on all best practices taught throughout the semester.
- Code organization and coding best practices will be graded.
- State diagram will be graded based on completeness, accuracy and ease of interpretation.
- Testing and analysis technical report will be graded based on presentation, completeness, and accuracy.

## How to Ask for Help

1. If you have a general / non-coding question, you should ask your TAs / Dr. Palmeri on Ed to allow any of them to respond in a timely manner.
2. Push your code to your GitLab repository, ideally with your active development on a non-main branch.
3. Create an [Issue](#) in your repository.
  - Add as much detail as possible as to your problem, and add links to specific lines / section of code when possible.
  - Assign the label “Bug” or “Question”, as appropriate.
  - Be sure to specify what branch you are working on.
  - Assign the Issue to one of the TAs.
  - If your TA cannot solve your Issue, they can escalate the Issue to Dr. Palmeri.
4. You will get a response to your Issue, and maybe a new branch of code will be pushed to help you with some example syntax that you can use `git diff` to visualize.

## What to Submit

- Make sure that all of your branches have been merged into the `main` branch.
- Create an annotated tag called `v1.0.0` to mark the commit that you want to be graded.
  - If you fix any bugs after creating this tag, you can create another tag called `v1.0.1`, etc.
  - Your latest tag will be the one that is graded up until the final due date/time of the project.
- Create an Issue in your repository with the title “Final Project Submission”, and assign it to Dr. Palmeri.
- **All repositories will be cloned at the due date/time for grading. Absolutely no changes will be accepted after this time.**

## Resources

### Heart Rate Service (GATT)

- [BLE Sample: Peripheral Heartrate](#)
- [Zephyr Docs: BT Heartrate Service](#)